

## **A protocol for rapid detection of the 2019 novel coronavirus SARS-CoV-2 using CRISPR diagnostics: SARS-CoV-2 DETECTR**

**v1: February 15, 2020**

**v2: February 18, 2020:** contact email provided, corrected RNase P gRNA sequence, updates to test interpretation matrix, corrections to LoD and assay reaction time for CDC SARS-CoV-2 qRT-PCR assay

**v3: March 2, 2020:** updates to minimum sample equipment and test interpretation matrix

James P. Broughton<sup>1</sup>, Wayne Deng<sup>2,3</sup>, Clare L. Fasching<sup>1</sup>, Jasmeet Singh<sup>1</sup>, Charles Y. Chiu<sup>2,3,4</sup>, Janice S. Chen<sup>1</sup>

<sup>1</sup>Mammoth Biosciences, Inc., South San Francisco, California, USA

<sup>2</sup>Department of Laboratory Medicine, University of California, San Francisco, California, USA

<sup>3</sup>UCSF-Abbott Viral Diagnostics and Discovery Center, San Francisco, California, USA

<sup>4</sup>Department of Medicine, Division of Infectious Diseases, University of California, San Francisco, California, USA

Acknowledgements: We thank Vikram Joshi, Nefeli Tsaloglou and Xin Miao for advice and helpful discussions in the preparation of this white paper.

Conflicts of Interest: JPB, CLF, JS and JSC are employees of Mammoth Biosciences, CYC is on the Scientific Advisory Board of Mammoth Biosciences, and JSC is a co-founder of Mammoth Biosciences. JPB, CLF, JS, CYC and JSC are co-inventors on CRISPR-related technologies.


Please contact [diagnostics@mammothbiosci.com](mailto:diagnostics@mammothbiosci.com) with any questions or feedback.

**\*\*\*DISCLAIMER: This protocol has not been approved by the FDA and should not be used as a clinical diagnostic\*\*\***


## Introduction

Given the global health emergency, rapid transmission, and severe respiratory disease associated with the outbreak of the 2019 novel coronavirus (SARS-CoV-2), Mammoth Biosciences has reconfigured our [DETECTR](#) platform to rapidly and accurately detect SARS-CoV-2 using a visual lateral flow strip format within 30 minutes from sample to result. To ensure specificity of detection, we selected a high-fidelity CRISPR detection enzyme and designed sets of gRNAs that can either 1) differentiate SARS-CoV-2 or 2) provide multi-coronavirus strain detection. SARS-CoV-2 DETECTR couples CRISPR detection with isothermal pre-amplification using primers based on protocols validated by the US Centers for Disease Control and Prevention (CDC) and World Health Organization (WHO). Currently in the United States, the CDC SARS-CoV-2 real-time RT-PCR diagnostic panel has a laboratory turnaround time of approximately 4-6 hours, with results that can be delayed for >24 hours after sample collection due to shipping requirements. In addition, these tests are only available in CDC-designated public health laboratories certified to perform high-complexity testing.

Mammoth is working to enable point of care testing (POCT) solutions that can be deployed in areas at greatest risk of transmitting SARS-CoV-2 infection, including airports, emergency departments, and local community hospitals, particularly in low-resource countries. Leveraging an “off-the-shelf” strategy to enable practical solutions within a short time frame, we describe here a protocol that is fast (<30 min), practical (available immediately from international suppliers), and validated using contrived samples.



**Table 1** | SARS-CoV-2 DETECTR assay workflow and specifications.



**Figure 1 | SARS-CoV-2 DETECTR has a limit of detection (n=7) of 156-625 copies per 20  $\mu$ l reaction (or 70-300 copies per  $\mu$ l input) and generates a clear visible signal on lateral flow strips within 30 minutes sample to result.**

## SARS-CoV-2 DETECTR Reagents

### Step 1: Isothermal amplification (62°C, 20 min)

[RT-LAMP Master Mix \(Supplier: NEB\)](#)


[DNA oligos \(Supplier: IDT\)](#)

Primer sequences:

Name	Sequence (5' → 3')
N-gene F3	AACACAAGCTTTCGGCAG
N-gene B3	GAAATTTGGATCTTTGTCATCC
N-gene FIP	TGCGGCAATGTTTGTAATCAGCCAAGGAAATTTGGGGAC
N-gene BIP	CGCATTGGCATGGAAGTCACTTTGATGGCACCTGTGTAG
N-gene LF	TTCCTTGTCTGATTAGTTC
N-gene LB	ACCTTCGGGAACGTGGTT
E-gene F3	CCGACGACGACTACTAGC



## Minimum sample equipment



- (1) DETECTR reagents
- (2) 37°C heat block
- (3) 62°C heat block
- (4) Nuclease-free water
- (5) Pipette tips
- (6) Pipette
- (7) Lateral flow strips

## Experimental Protocol

### A. Prepare nucleic acid sample and CRISPR reagents

1. Extract patient RNA following [CDC recommendations](#).
2. Prepare LbCas12a RNP complexes for the samples to be tested. One complex for N-gene, E-gene, and RNase P gRNAs is needed for each sample.

Reagent	Volume	Final Concentration
Nuclease-free water	15.75 $\mu$ L	
10X NEBuffer 2.1	2 $\mu$ L	1X
1 $\mu$ M LbCas12a	1 $\mu$ L	50 nM
1 $\mu$ M gRNA	1.25 $\mu$ L	62.5 nM
<b>TOTAL VOLUME</b>	<b>20 <math>\mu</math>L</b>	

- Incubate LbCas12a with gRNA to generate RNP complexes for 30 minutes at 37°C.
- Add reporter substrate to final concentration of 500 nM.
- Place reactions on ice until ready to proceed.
  - Complexes are stable at 4°C for at least 24 hours.

#### B. Run DETECTR reaction

- On ice, prepare three RT-LAMP reactions, one each for N-gene, E-gene, and RNase P primer sets:

Reagent	Volume	Final Concentration
10X Isothermal Amplification Buffer (NEB)	2.5 $\mu$ L	
100 mM MgSO <sub>4</sub> (NEB)	1.13 $\mu$ L	6.5 mM (4.5 mM added, 2 mM in 1X IsoAmp Buffer)
10 mM dNTPs (NEB)	3.5 $\mu$ L	1.4 mM
10X Primer Mix	2.5 $\mu$ L	0.2 $\mu$ M F3 0.2 $\mu$ M B3 1.6 $\mu$ M FIP 1.6 $\mu$ M BIP 0.8 $\mu$ M LF 0.8 $\mu$ M LB
Bst 2.0 polymerase (NEB)	1 $\mu$ L	8 units / rxn
Warmstart RTx (NEB)	0.5 $\mu$ L	7.5 units / rxn
Nuclease-free water	3.87 $\mu$ L	

Nucleic acid sample	5 $\mu$ L	
<b>TOTAL VOLUME</b>	<b>25 <math>\mu</math>L</b>	

2. Incubate at 62°C for 20 minutes.
  - a. Note: Use precaution when opening amplification tubes to prevent amplicon contamination.
3. Combine 2  $\mu$ L of the RT-LAMP reaction with 20  $\mu$ L of the LbCas12a RNP complex with the appropriate gRNA.
4. Add 80  $\mu$ L of 1X NEBuffer 2.1.
5. Incubate at 37°C for 10 minutes.
6. Insert Milenia HybriDetect 1 (TwistDx) lateral flow strip directly into reaction.
7. Allow the lateral flow strip to run for 2 minutes at room temperature and observe the result.


### C. Test interpretation

Note: The line closest to the sample pad is the control line and the line that appears farthest from the sample pad is the test line (see Figure 1). A sample with complete cleavage of the reporter molecule may appear to have no signal at the control line.

N-gene	E-gene	RNase P	Interpretation
+	+	+/-	SARS-CoV-2 positive
+	-	+/-	Indeterminate
-	+	+/-	Indeterminate
-	-	+	SARS-CoV-2 negative
-	-	-	QC failure

## Appendix

While we were preparing this white paper, another [protocol for SARS-CoV-2 detection using CRISPR diagnostics \(SHERLOCK, v.20200214\)](#) was published. We compare the assay workflows and specifications between CRISPR diagnostics and established CDC/WHO protocols below. (Note: as of this publication, CRISPR diagnostics workflows have not yet been approved by the FDA)



**Appendix Figure 1** | Comparison of SARS-CoV-2 assay workflows for DETECTR, SHERLOCK, and CDC/WHO.

	<b>SARS-CoV-2 DETECTR</b>	<b>SARS-CoV-2 SHERLOCK</b>	<b>CDC SARS-CoV2 qRT-PCR</b>
Target	N gene & E gene  (N gene gRNA compatible with CDC N2 amplicon, E gene compatible with WHO protocol)	S gene & Orf1ab gene	N-gene (3 amplicons)
Sample control	RNase P	None	RNase P
Limit of Detection	70-300 copies/μl input	10-100 copies/μl input	3.16-10 copies/μL input
Assay reaction time	~30 min	~60 min	~120 minutes



Assay components	RT-LAMP (62°C, 20 min) Cas12 (37°C, 10 min) Lateral flow (RT, 2 min)	RT-RPA (42°C, 25 min) IVT + Cas13 (37°C, 30 min) Lateral flow (RT, 2 min)	UDG digestion (25°C, 2 min), reverse transcription (50°C, 15 min), denature (95°C, 2 min) amplification, (95°C, 3 sec; 55°C 30 sec; 45 cycles)
Heavy instrumentation required	No	No	Yes
FDA EUA approval	No	No	Yes

**Appendix Table 1 |** Comparison of SARS-CoV-2 specifications for CRISPR diagnostic protocols to the current CDC assay.